AUSTRIA

Bauerfeind Ges.m.b.H.
Hainburger Straße 33
1030 Wien
P +43 (0) 800 4430-130
F +43 (0) 800 4430-131
E info@bauerfeind.at

BENELUX

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 english

Dear Customer,

thank you for choosing a Bauerfeind product.

We work to improve the medical effectiveness of our products every day – because your health is very important to us. Please closely observe the specifications in these instructions for use and the instructions given by the specialist. If you have any questions, please contact your doctor or medical retailer.

Intended purpose

SofTec Genu is a medical device. It is an orthosis¹ for active and passive stabilization of the knee joint – ideal for long-term treatment.

Indications

- Rupture of the anterior and / or posterior cruciate ligament
 - Non-surgical treatment
 - Pre-operative and delayed postoperative management
 - Long-term rehabilitation
 - Sports therapy
 - Chronic insufficiency
- Severe and / or complex instability (traumatic, degenerative, post-operative, e.g. "unhappy triad", endoprosthetics, Genu recurvatum)
- Non-surgical functional treatment of collateral ligament lesions
- Treatment of an "unstable knee joint"
 - Severe osteoarthritis of the knee
 - Rheumatoid arthritis (PCP)
- To restrict the range of movement of the knee (e.g. after meniscal suturing or meniscal implantation)

Risks of using this product

Important information

- The SofTec Genu must be individually adjusted for the patient by a trained specialist² in order to ensure maximum effectiveness and optimum wearing comfort¹.
- SofTec Genu is a prescribed product that should be applied under a physician's guidance. SofTec Genu should only be applied in accordance with the instructions contained in these usage directions and for the listed areas of application.
- Before using the product, the specialist or your physician must determine whether thrombosis prevention is necessary.

- Before applying SofTec Genu after an operation, be sure that the sensitivity and swelling are reduced enough to ensure that there are no unnecessary pressure points or complications. Typically, SofTec Genu is not applied post-operatively for 5–8 days, depending on the individual patient's pain and swelling.
- Modifications of SofTec or any joint bar adjustments may only be made by suitably trained specialists². Failure to comply with this requirement may adversely affect product performance, thereby cancelling any product liability.
- If SofTec Genu is used during exercise or sporting activities do not limit the flexion of the hinge. It is possible to limit extension to 0°, to avoid overstretching the knee joint.
- SofTec Genu is not designed to prevent injuries, especially when engaging contact sports. Contact sports inherently pose a greater risk of injury.
- A combination with other products (e.g. compression stockings) must be discussed beforehand with your medical advisor.
- When resting, the straps can be loosened for better comfort.
- Avoid the use of ointments, lotions or any substances including grease or acids.
- Side effects, that affect the entire organism, have not been reported to date. Correct fitting is assumed. Any supports and orthoses applied externally to the body can, if tightened excessively, lead to local pressure signs or, rarely, constrict the underlying blood vessels or nerves. In such cases, please loosen the straps slightly and have the correct fit and size of the orthosis checked if necessary.
- If you notice any changes or an increase in symptoms while wearing the product, stop any further use and contact your physician.

 Contains natural rubber latex (elastodiene), which can cause skin irritation in persons with an allergy. Take off the orthose if you notice any allergic reaction.

Contra-indications

Hypersensitive reactions harmful to health have not been reported to date. In the following conditions such aids should only be fitted and worn after consultation with your medical advisor:

- Skin disorders / injuries in the relevant part of the body, particularly if inflammation is present. Likewise, any raised scars with swelling, redness and excessive heat build-up
- Venous insufficiency, particularly with varicose veins
- Impaired sensation and circulatory disorders of the legs / feet, e.g. diabetes mellitus
- Impaired lymph drainage – including soft tissue swellings of uncertain origin located away from the fitted aid
- Post-thrombotic syndrome

Application instructions

Application and Removal of SofTec Genu:

- Open the zipper and all the straps. Pull the straps out of the plastic loops. Fasten down the ends of the straps on the Velcro area.
- Hold SofTec Genu with both hands at the upper end (with the label). Turn SofTec Genu so that the rear part (zipper) points to the front. (Fig. 1).
- Pull SofTec Genu over your ankle until the lower end is above your ankle.
- Turn SofTec Genu by 180° (half a turn) around your leg. (Fig. 2).
- Pull SofTec Genu up your leg until your kneecap is centered in the shaped ring (Fig. 3).
- The zipper can be closed easily by pulling the front lower part of SofTec Genu in the area of the zipper to the back until it overlaps. (Fig. 4).
- Guide all the straps through the plastic loops and fasten them loosely. An intended pressure is caused by tightening the straps in the following order. (Fig. 5):
 - Calf strap furthest from the knee (1)
 - Thigh strap closest to the knee (2)
 - Thigh strap furthest from the knee (3)
 - Calf strap closest to the knee (4).
- Be sure SofTec Genu is seated properly. Any straps that are too tight should be loosened. If you get pressure sores in the area of the hinges on the sides, get your medical retailer to reshape them.
- Remove SofTec Genu by doing the above strapping sequence in the reverse order.

Cleaning instructions

Never expose the orthosis to direct heat (e.g. heater, sunlight, etc.)!

SofTec Genu can be washed separately by hand with a mild detergent at a temperature of 30 °C (86 °F). The hinged side bars must be removed from the side pockets. This is done as follows:

- Open the two thigh straps. Pull the straps out of the plastic loops. Fasten down the ends of the straps on the Velcro area. The side pockets can be opened by pulling the entire Velcro section to one side. (Fig. 6).
- Remove the hinged side bars.  **Important:** bear in mind that the joint bars are contoured to fit the body and are therefore different. The hinged side bars and the corresponding side pockets are marked with the same color so that the hinged side bars can be replaced in the proper side pocket.
- Before washing, guide all the straps through the plastic loops and close the Velcro fastenings in order to preserve their adhesive strength over a long period of time.
- After SofTec Genu is washed and air dried, replace the hinged side bars into the proper side pocket. Be sure the colors on the hinged side bars and side pockets match. The upper side bar must be placed in the inner pocket. (Fig. 7).
- Close the side pockets. (Fig. 8).
- After strenuous activity, we recommend removing the joint splints, wiping them with a damp cloth and leaving them to dry separately.

This product has been examined by our own quality control system. However, should you have any complaint, please contact your retailer or our customer service department.

Part of the body this product is used for

In accordance with the indications (knee). See purpose.

Maintenance instructions

If handled and cared for correctly, the product is practically maintenance-free.

Technical data

The SofTec Genu is an orthosis for the knee joint. It consists of two lateral hinges, four straps, knitted fabric and a pad.

Notes on reuse

This product is intended to be fitted to one patient only.

Warranty

The statutory regulations of the country of purchase apply. Please first contact the retailer from whom you obtained the product directly in the event of a potential claim under the warranty. The product must be cleaned before submitting it in case of warranty claims. Warranty may be limited or excluded if the instructions on how to handle and care for the SofTec Genu have not been observed.

Warranty is excluded if:

- The product was not used according to the indication
- The instructions given by the specialist were not observed
- The product was modified arbitrarily

Disclaimer

Do not diagnose yourself or decide on which medication to take unless you are a medical specialist. Before first using our medical product, please seek the advice of a physician or trained specialist, as this is the only way to assess the effect of our product on your body and determine the potential risks of using this product due to your personal constitution. Follow the advice of this specialist as well as all instructions in this document or its online version, including any excerpts (as well as texts, images, graphics etc.). If you still have doubts after your consultation with the specialist, please contact your physician or retailer, or contact us directly.

Duty to report

Due to regional legal regulations, you are required to immediately report any serious incident involving the use of this medical device to both the manufacturer and the responsible authority. Our contact details can be found on the back of this brochure.

Disposal

Once you have finished using the product, please dispose of it according to local specifications.

Material content

Polycarbonate (PC), Polyamide (PA), Silicone (SI), Polyurethane (PUR), Spandex (EL), Polyoxymethylene (POM), Elastodiène (ED), Polyester (PES), Cotton (CO), Zinc alloy, Polyethylene, high density (PE-HD)

MDI – Medical Device

UDI – DataMatrix Unique Device Identifier UDI

Version: 2022-09

¹Orthosis = orthopaedic appliance used to stabilize, relieve, immobilize, control or correct the limbs or the torso

²A specialist is any person who is authorized according to the state regulations for fitting and instruction in the use of supports and orthoses which are relevant to you.

Specialists²

Application instructions

Size calculation

Please determine the correct size in accordance with the specifications on the packaging.

Assembly and fitting instructions

Procedure for First Time Application Contouring the Hinged Side Bars

- Open the two thigh straps. Pull the straps out of the plastic loops and fasten down the ends of the straps on the Velcro area. You can now open the side pockets and remove the hinged side bars (Fig.1).
- Check the fit of the joint bars on the inner and outer aspects of the knee and reshape if necessary.

Important!

- When contouring the sidebars be sure not to damage the hinges. The use of bending irons, with protection if necessary, is recommended for preventing such damage. To correctly fit the brace, undo the adjusting screws (Fig.10), making sure that the stop device is not lifted out of its guide. Avoid bending the sidebars too close to the hinge. The hinged side bars must be contoured so that no pressure is exerted at the level of the top of the calfbone (top of the fibula) to prevent irritation.
- Put the hinged side bars back into the corresponding side pockets. The joint bars for the inner (medial) and outer (lateral) aspects of the knee are marked with different colors. The color markings on the side bar and the color of the inner side pockets must match. The side bars must not be pushed under the colored inner pockets at this point (Fig.2).
- Close the side pocket.

First Time Application of the SofTec Genu

- Open the zipper and all the straps. Fasten down the ends of the straps on the Velcro area.
- Hold SofTec Genu with both hands at the proximal end (near the label). Turn SofTec Genu so that you can see the zipper (Fig. 3).
- Pull SofTec Genu over the ankle until the distal end is above the ankle.
- Turn SofTec Genu by 180° (half a turn) around the leg (Fig.4).
- Pull the SofTec Genu over the leg until the kneecap is in the center of the pad (Fig.5).
- The zipper can be closed easily by pulling the front lower part of SofTec Genu in the area of the zipper to the back until it overlaps (Fig.6).
- Guide all the straps through the plastic loops and fasten them loosely. The intended pressure is achieved by tightening the straps in the following order (Fig.7):
 - Distal calf strap (1)
 - Proximal thigh strap (2)
 - Distal thigh strap (3)
 - Proximal calf strap (4)
- Be sure SofTec Genu is seated properly. Any straps that are too tight should be loosened. If the patient feels any pressure points in the area of the hinged side bars they may need to be contoured again.

Setting the Patients Axis of Rotation

The hinges of SofTec Genu come with an initial setting of 16 mm rearward offset. This is obvious if the unit is set for the middle (third) embossed mark to be parallel to the proximal (upper) area of the distal (lower) side bars when the side bar is bent by 30° (Fig.8). The hinges have a technically intelligent mechanism that allows optimal adjustment to the patient’s axis of rotation. The individual patients axis of rotation can be determined as follows:

- Open the distal (lower) thigh strap, fasten down the end, and open the lower area of the side pockets (Fig.9).
- Find the stop screws that are located on the proximal (upper) side bars, just above the hinge. Unscrew them until they stop.
- Press downward on the stop screws and move them in the direction of the knee cap (anteriorly) until they stop (Fig.10). (To check only slightly rotate the plastic covered portion.) Check that the default setting of 16 mm rearward offset still applies after opening the stop. This can be seen from the positioning of the metal pins in the locator (Figs.11 and 12). Tighten the two stop screws.
- Close the distal (lower) thigh strap (loosely). Have the patient flex and extend the knee joint in a range of motion between 10° and 45°, approximately 15 times. This is necessary to find the individual patient’s axis of rotation. If the patient is not able to bear weight after the operation, it is only necessary to swing the patient’s knee in a range of motion between 10° and 45°. As the patient becomes more active, it is necessary to readjust the axis of rotation under weight bearing conditions.
- Open the distal thigh strap, fasten down the end and open the lower area of the side pockets (Fig.13).
- Unscrew the stop screws and push them towards the back of the kneecap (posteriorly) until they click into place (Fig.14). If they do not immediately click into place, slightly rotate the plastic covered portion of the hinge (max. 5°) (Fig.15).
- Tighten both stop screws again. Now the patient’s individual axis of rotation has been found and secured.
- Open the proximal thigh strap and the side pockets.
- Push the upper portion of the side bars under the colored inner pockets (Fig.16).
- Close the two side pockets and the thigh straps.
- To ensure the intended pressure, close the straps in the following order (Fig.17):
 - Distal thigh strap (1)
 - Proximal thigh strap (2)

Important

The individual patients axis of rotation should be readjusted if the circumference of the patient’s leg changes in the course of therapy, e.g., as a result of muscle build-up. In such a case you should repeat the procedure described above from the section “Finding the patients axis of rotation”.

Setting the range of motion

The flexion and extension of the joints on SofTec Genu can be restricted in 10° increments.

- Open the two thigh straps and fasten the ends down. The orthosis pocket on the outer aspect of the knee can now be opened down to the joint mechanism (Fig.18).
- Remove the joint cover on the outer aspect of the knee.
- Bend the lower leg until the embossed marks on the proximal (upper) side bar become parallel to the upper portion of the distal (lower) side bar (30°) (Fig.19).
- You will find 2 metal pins with ring-type handles located in holes in the middle of the hinge. Pull them out.
- Rotate the scaling disc until the arrow on the upper portion of the distal side bar lines up with the arrow on the scaling disc (Fig.20). Select the flexion and extension angles and place the metal pins in the corresponding hole that most correctly lines up with the numbers. Each hole is a 10° increment.
- Replace the plastic hinge cover so that the pins stay in the holes.
- Now open the orthosis pocket on the inner aspect of the knee.
- Remove the plastic hinge cover.
- Remove the metal pins by the ring-type handles.
- Flex the knee gently until the hinge stops. Place a pin in the appropriate hole; it does not necessarily have to match with the lateral side. Now extend the knee until the hinge stops and place the pin in the hole that matches the lateral side. Important! The settings for the restriction of flexion and extension on the inner and outer aspects of the knee must be identical.

- Replace the plastic hinge cover and close the medial side pocket and the thigh straps.
- To secure SofTec Genu properly the strapping sequence should be as follows (Fig.21):
 - Distal thigh strap (1)
 - Proximal thigh strap (2)

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